

Replacement Recline Sling

F15

- Designed to provide closed leg full body support on transfer
- Gives reclined patient transfer position
- Pommel strap for additional comfort and safety
- To be applied in prone position
- Soft binding as standard
- Fitted colour coded Silva Safety Slot tapes
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 71°C/159°F
- Tumble dry low heat (max 50°C)



Product Codes & Specifications:



Silva-Superfine PLUS®		A	B	C	D
XXS	F15SFPXXS	66 - 76cm	78 - 81cm	60cm	69cm
XS	F15SFPXS	76 - 86cm	81 - 92cm	70cm	75cm
S	F15SFPS	86 - 96cm	92 - 104cm	78cm	83cm
M	F15SFPM	96 - 108cm	104 - 116cm	88cm	93cm
L	F15SFPL	108 - 130cm	116 - 128cm	108cm	106cm
XL	F15SFPXL	130 - 144cm	128 - 140cm	118cm	117cm
XXL	F15SFPXXL	144 - 160cm	128 - 140cm	132cm	137cm

