





Tinkham Sling with HS

- Fully supportive sling suitable for patients who have difficulty in maintaining body control
- Divided leg sling
- Soft padded neck roll and head support
- Colour coded, Silva-Safety Slot® shoulder and leg tapes supplied as standard
- Two sets of adjustable, sliding hip tapes supplied as standard
- Features a wrap-around thoracic section, designed to fit firmly around the body, with clunk-click buckle
- Wash temp 71°C / 159°F
- Tumble dry low heat (max 50°C)

Product Codes & Specifications:



	Silva-Superfine PLUS®	A	B	C	D
XXXS	TINKSFPHS-XXXS	56 - 66cm	50 - 55cm	35cm	51cm
XXS	TINKSFPHSXXS	66 - 76cm	56 - 62cm	39cm	61cm
XS	TINKSFPHSXS	76 - 86cm	60 - 66cm	49cm	71cm
S	TINKSFPHSS	86 - 96cm	66 - 74cm	57cm	77cm
M	TINKSFPHSM	96 - 108cm	74 - 80cm	69cm	89cm
L	TINKSFPHSL	108 - 130cm	80 - 88cm	77cm	102cm
XL	TINKSFPHSXL	130 - 144cm	88 - 100cm	85cm	108cm



TINKSFPHS



NB: Image shows colour coded tapes, this product is supplied with colour coded Silva-Safety Slot® tapes as standard.