


CATALYST 5 & 5T





INNOVATIVE AND PATENTED FRAME DESIGN

Patented frame concept uses intelligent geometry to create a significantly stronger, lighter and better performing chair.

STANDARD & CURVED AXLE PLATES

Catalyst 5 provides the ability to adjust CG over a 3 1/2" range. The adjustment is a snap with the innovative 3D locking index washer design and requires only one wrench.

SWING AWAY FOOTRESTS

Enjoy the simplicity and strength of Catalyst's patented swing-away footrest hangers. The tapered engagement keeps the footrests snug and rattle free. Comes with a lifetime warranty.

UNMATCHED PERFORMANCE.

The professional's best choice in folding ultra-lightweight mobility.

Also Available in Titanium!



SPECIFICATIONS

5 Standard Weight:	23 lbs.*	10.43 kg*
5Ti Standard Weight:	21 lbs.*	9.53 kg*
Seat Width:	14-22 in.	36-56 cm
Seat Depth:	14-20 in.	36-51 cm
Front Seat Height:	13.5-21 in.	34-53 cm
Rear Seat Height:	11.5-19.5 in.	29-50 cm
Top Frame Tubing Size:	1 in.	2.54 cm
Bottom Frame Tubing Size:	1.125 in.	2.86 cm
Weight Capacity:	250/350 lbs.	113/158 kg
Ansi Resna WC-19 Transit:	Yes	Yes

* With spoke wheels and swing-away armrests; without footrests.