

Little Wave "Click"



Not Just a Wheelchair. A Member of the Team.

DESIGNED FOR OUR *NEXT* GENERATION.



Welcome to the most intuitive pediatric wheelchair you've ever seen.

An adjustment system designed to provide accurate, easy growth options. A system so easy to adjust it promotes proper positioning and wheel access so your child gains the active independence and mobility they want.

For the first time ever, growth is just a "Clik" away.

The Clik's proprietary **Index System** easily indexes and aligns as you optimize your adjustment. Providing you just enough guidance to achieve the very best results, every time. And the confidence to know it's right.





Evolving and Growing with Your Child.

Seat height, width, depth, CG, caster wing position, seat upholstery, seat pans...the most built in adjustment available. With all the possibilities, meeting kids needs has never been so easy.

YOU ARE ONLY A "CLIK" AWAY,
FROM JUST RIGHT.





Expanding your child's world with confidence.

BE DYNAMIC.

Being able to move yourself in space is an essential phase of neurodevelopment. The challenge is maintaining security and stability. Our new adjustable dynamic 5th wheel, that allows you to adjust both the range and spring rate, makes it possible.

Whether setup in mid-wheel for the smallest riders or to help promote and teach the balance and wheeling skills that encourage confidence to overcome everyday obstacles and terrain. **Be Dynamic.**



Kliq's unique swing away armrest allows for independent angle adjustability, while integrated into the backrest for easy growth. Your arm is always where you need it.




The replaceable cross tube with TAPER LOK increases Kliq's rigidity resulting in an easier to push chair. Optional integrated receivers use less space so you can get everything in the right spot. The first time.



Compound offset frames ensure your kid has the fit they need. More room for abducted positioning, AFO's and braces.



WHEELS, SPOKE GUARDS, FRAMES & ACCENTS.
THOUSANDS OF FUN COMBINATIONS!



SPECIFICATIONS

Transport Weight:	12.5 lbs.	5.67 kg
Seat Width:	8-16 in.	20-41 cm
Seat Depth:	8-16 in.	20-41 cm
Front Seat Height:	11.5-20 in.	29-51 cm
Rear Seat Height:	11.5-20.5 in.	29-52 cm
Tubing Size:	1.125 in.	2.85 cm
Weight Capacity:	165 lbs.	74.8 kg
WC-19 Transit:	Yes	Yes

