

ROGUE



Ki Mobility
Better by Design



DESIGN

You will not see anything else like it. Minimalist in everything but the frame tubing, Rogue catches the eye. Incredibly light at 10.6 lbs. Rogue conserves your energy. An integrated locking back mechanism makes transportation easier than it has ever been before.

RIDE

You will not feel anything else like it. Strategically designed with oversized 1 3/8" 7000 series tubing the Rogue has the ideal combination of vibration dampening and rigidity. It's the most responsive rigid frame ride on the market.

PRECISION

It is the only indexed caster adjustment system that gives a perfect alignment for everyone at any seat angle. The modified rack and pinion gear is infinitely adjustable over any range and can be adjusted while the chair is occupied. There are no small components to take apart and put back together. It is as simple as tightening and loosening a screw. It is called RP Tuning. And only Rogue has it.

It's in the soul of the Rogue, the way it stands out, its defiance of convention, its independence that separates it from all others.

RESIST CONFORMITY. GO ROGUE.



SPECIFICATIONS

| | | |
|---------------------------|---------------|----------|
| Transport Weight: | 10.6 lbs.* | 4.81 kg* |
| Seat Width: | 12-20 in. | 30-51 cm |
| Seat Depth: | 14-20 in. | 36-51 cm |
| Front Seat Height: | 13.5-20.5 in. | 34-52 cm |
| Rear Seat Height: | 12.5-20 in. | 32-51 cm |
| Tubing Size: | 1 3/8 in. | 3.49 cm |
| Weight Capacity: | 275 lbs. | 124 kg |
| Ansi Resna WC-19 Transit: | Yes | Yes |

* Rogue weight calculation based on: 16"x16" Standard Length Frame, Performance Sling Position, Adjustable Height Back - 12", 4" Poly Casters, Tubular Open Footplate and Standard Upholstery.