

Please mark clearly with an X the items you wish to purchase.
 Email form to NuVision Rehab Group at sales@nuvisionrehab.com or fax form to 403-253-1744.
 Contact our customer service team at (877) 252-6655.

Account #: _____	Ship to: _____
Date: _____	Name: _____
PO #: _____	Address: _____
Buyer: _____	City: _____
Attention: _____	Phone: () _____

EIO Push Chair Features:

- Removable hip lateral and trunk lateral supports
- Height adjustable padded headrest
- Back rest angle adjusts from 25° to 75°
- Footwell and angle adjustable push handle
- Adjustable 5-Point padded chest harness and shoulder pads as well as a padded 3-point locking pelvic belt
- Removable sun canopy with window and pocket
- Solid filled 9" front and 11" rear wheels require no air
- Swivel front wheels with lock-out feature add security on rough terrain
- Angle adjustable foot well with a 2-position foot-block.
- Rear spring suspension for comfortable ride
- Easy access one-touch parking brake
- Lightweight aluminum frame
- Shopping basket included

2 year warranty
 EIO Push Chair

5 year warranty
 Liners/Sitters



Includes
 2-Position
 Foot-Block

Optional Accessories:



Rain Cover



For kids who need mild to moderate positioning, the Size 1 and 2 Soft-Touch® Liners and Size 1 and 2 Soft-Touch® Sitters attach to the stroller.



EIO Push Chair







Please mark clearly with an X the items you wish to purchase.
 Email form to NuVision Rehab Group at sales@nuvisionrehab.com or fax form to 403-253-1744.
 Contact our customer service team at (877) 252-6655.





Push Chair Specifications	Stroller Only	Liner Size 1	Liner Size 2	Sitter Size 1	Sitter Size 2
Seat to Footwell	13"	15"	15"	16"	17"
Seat	14"W x 12"D				
Back Height	24"				
Maximum User Weight	90 lbs.				
Stroller Weight	22 lbs.				
Handle Height to Floor	30" - 43"				
Dimensions when Open	45"H x 24"W x 46"D				
Dimensions when Folded	16"H x 18"W x 32"D				

Item	Part #	MSRP
Special Tomato® EIO Push Chair	_____ 63000000	\$ 1,455.23
Rain Cover	_____ 63030000	\$ 119.28




Soft-Touch® Liners Back Cushions

Size	Cherry 	Lilac 	Aqua 	Dark Gray 	MSRP
LB1	_____ 72050101	_____ 72050102	_____ 72050104	_____ 72050109	\$ 315.64
LB2	_____ 72050201	_____ 72050202	_____ 72050204	_____ 72050209	\$ 343.06

Soft-Touch® Liners Seat Cushions

Size	Cherry 	Lilac 	Aqua 	Dark Gray 	MSRP
LS1	_____ 72040101	_____ 72040102	_____ 72040104	_____ 72040109	\$ 315.64
LS2	_____ 72040201	_____ 72040202	_____ 72040204	_____ 72040209	\$ 343.06

Soft-Touch® Sitters

Size	Lilac 	Aqua 	Dark Gray 	MSRP
Size 1	_____ 77010102	_____ 77010104	_____ 77010109	\$ 660.10
Size 2	_____ 77010202	_____ 77010204	_____ 77010209	\$ 868.58

Toll-Free: (877) 252-6655 Fax: (403) 253-1744 Email: info@nuvisionrehab.com
 Order Forms Available Online at: www.nuvisionrehab.com