

Please mark clearly with an X the items you wish to purchase.
 Email form to NuVision Rehab Group at sales@nuvisionrehab.com or fax form to 403-253-1744.
 Contact our customer service team at (877) 252-6655.


Account #: _____	Ship to: _____
Date: _____	Name: _____
PO #: _____	Address: _____
Buyer: _____	City: _____
Attention: _____	Phone: () _____

Height Right™ Chair Features:

- Adjustable footrest for proper positioning
- Seat and footrest adjust allowing the chair to grow with the child
- Optional Soft-Touch® cushions easily attach for children who need lateral support
- Optional Padded Seat and Back Cushions are designed to fit the Height Right™ Chair.



Optional Accessories:





Height Right™ Chair







Please mark clearly with an X the items you wish to purchase.
 Email form to NuVision Rehab Group at sales@nuvisionrehab.com or fax form to 403-253-1744.
 Contact our customer service team at (877) 252-6655.

Height Right™ Chair Specifications	
Seat Plate	16.5"W × 10"D
Seat Plate Height Range (from Floor)	16.5" - 24"
Foot Plate	11.5"W × 16.5"D
Foot Plate Height Range (from Floor)	4" - 15"
Tray	19"W × 8"D
Tray Height from Floor	26"
Overall Height	33"
Overall Width	17"
Overall Depth	21"
Weight Capacity	250 lbs.

Item	Part #	MSRP
Height Right™ Chair	___ 61000000	\$ 404.44

Optional Accessories	Part #	MSRP
Wooden Tray (w/Protective Cover)	___ 61010000	\$ 167.77
Plastic Tray Cover ONLY	___ 61020000	\$ 96.16
Footcups (pair) - Dark Gray	___ 61060009	\$ 302.65

Optional Cushions	Cherry 	Lilac 	Aqua 	Dark Gray 	MSRP
Soft-Touch® Seat Cushion (HS)	___ 61040001	___ 61040002	___ 61040004	___ 61040009	\$ 302.65
Soft-Touch® Back Cushion (HB)	___ 61050001	___ 61050002	___ 61050004	___ 61050009	\$ 302.65

WARNING: This chair is as stable as a typical kitchen chair or baby highchair. However, if your child self-stimulates using vestibular movements (rocking back and forth), this product might not be the best choice to meet the child's needs.