

Please mark clearly with an X the items you wish to purchase.  
Email form to NuVision Rehab Group at [sales@nuvisionrehab.com](mailto:sales@nuvisionrehab.com) or fax form to 403-253-1744.  
Contact our customer service team at (877) 252-6655.

Account #: \_\_\_\_\_ Ship to: \_\_\_\_\_  
Date: \_\_\_\_\_ Name: \_\_\_\_\_  
PO #: \_\_\_\_\_ Address: \_\_\_\_\_  
Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
Attention: \_\_\_\_\_ Phone: (     ) \_\_\_\_\_

### Jogger Features:

- Super-light aluminium frame folds easily
- Comfortable padded seat reclines
- Adjustable 5-point padded chest harness and shoulder pads as well as a padded 3-point locking pelvic belt
- 12" pneumatic wheels with front swivel lock-out
- Removable machine-washable fabric

2 year warranty

Jogger

5 year warranty

Liners/Sitters



One-step, compact folding design  
(Canopy folds with Jogger)

Front wheel rotates 360° and includes lock-out pin for rough terrain



### Optional Accessories:



Optional rain cover



Soft-Touch liners and padded headrest



Soft-Touch Sitter



Removable canopy with window and height adjustable handle




Please mark clearly with an X the items you wish to purchase.  
 Email form to NuVision Rehab Group at sales@nuvisionrehab.com or fax form to 403-253-1744.  
 Contact our customer service team at (877) 252-6655.

Jogger Specifications	Stroller Only	Liner Size 1	Liner Size 2	Sitter Size 1	Sitter Size 2
Seat to Footwell	14"	16"	16"	17"	18"
Seat	14"W x 12"D				
Back Height	25"				
Maximum User Weight	110 lbs.				
Stroller Weight	27 lbs.				
Handle Height to Floor	35" - 43"				
Recline Range	100° - 140°				
Dimensions when Open	43"H x 26"W x 42"D				
Dimensions when Folded	16"H x 21"W x 33"D				





Item	Part #	MSRP
Jogger	___ 62000000	\$ 1,488.74

Optional Accessories	Part #	MSRP
Medical Necessities Bag (10 lb. weight capacity)	___ 62010000	\$ 119.28
Head Support	___ 62020000	\$ 227.19
Rain Cover	___ 62030000	\$ 119.28




#### Soft-Touch Liners Back Cushions

Size	Cherry 	Lilac 	Aqua 	Dark Gray 	MSRP
LB1	___ 72050101	___ 72050102	___ 72050104	___ 72050109	\$ 315.64
LB2	___ 72050201	___ 72050202	___ 72050204	___ 72050209	\$ 343.06

#### Soft-Touch Liners Seat Cushions

Size	Cherry 	Lilac 	Aqua 	Dark Gray 	MSRP
LS1	___ 72040101	___ 72040102	___ 72040104	___ 72040109	\$ 315.64
LS2	___ 72040201	___ 72040202	___ 72040204	___ 72040209	\$ 343.06

#### Soft-Touch Sitters

Size	Lilac 	Aqua 	Dark Gray 	MSRP
Size 1	___ 77010102	___ 77010104	___ 77010109	\$ 660.10
Size 2	___ 77010202	___ 77010204	___ 77010209	\$ 868.58