

Please mark clearly with an X the items you wish to purchase.  
Email form to Ki Mobility Canada Ltd. at [sales@kimobilitycanada.ca](mailto:sales@kimobilitycanada.ca) or fax form to 403-253-1744. Contact our customer service team at (877) 252-6655.

Account #: \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO #: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (     ) \_\_\_\_\_

### Jogger Features:

- Super-light aluminium frame folds easily
- Comfortable padded seat reclines
- Adjustable 5-point padded chest harness and shoulder pads as well as a padded 3-point locking pelvic belt
- 12" pneumatic wheels with front swivel lock-out
- Removable machine-washable fabric

**2** year warranty

Jogger

**5** year warranty

Liners/Sitters



One-step, compact folding design  
(Canopy folds with Jogger)

Front wheel rotates 360° and includes lock-out pin for rough terrain



### Optional Accessories:



Optional rain cover



Soft-Touch liners and padded headrest



Soft-Touch Sitter



Removable canopy with window and height adjustable handle




Please mark clearly with an X the items you wish to purchase.  
 Email form to Ki Mobility Canada Ltd. at sales@kimobilitycanada.ca or fax form to  
 403-253-1744. Contact our customer service team at (877) 252-6655.

Jogger Specifications	Stroller Only	Liner Size 1	Liner Size 2	Sitter Size 1	Sitter Size 2
Seat to Footwell	14"	16"	16"	17"	18"
Seat	14"W x 12"D				
Back Height	25"				
Maximum User Weight	110 lbs.				
Stroller Weight	27 lbs.				
Handle Height to Floor	35" - 43"				
Recline Range	100° - 140°				
Dimensions when Open	43"H x 26"W x 42"D				
Dimensions when Folded	16"H x 21"W x 33"D				





Item	Part #	MSRP
Jogger	___ 62000000	\$ 1,849.01

Optional Accessories	Part #	MSRP
Medical Necessities Bag (10 lb. weight capacity)	___ 62010000	\$ 133.83
Head Support	___ 62020000	\$ 254.91
Rain Cover	___ 62030000	\$ 133.83




#### Soft-Touch Liners Back Cushions

Size	Cherry 	Lilac 	Aqua 	Dark Gray 	MSRP
LB1	___ 72050101	___ 72050102	___ 72050104	___ 72050109	\$ 354.14
LB2	___ 72050201	___ 72050202	___ 72050204	___ 72050209	\$ 384.91

#### Soft-Touch Liners Seat Cushions

Size	Cherry 	Lilac 	Aqua 	Dark Gray 	MSRP
LS1	___ 72040101	___ 72040102	___ 72040104	___ 72040109	\$ 354.14
LS2	___ 72040201	___ 72040202	___ 72040204	___ 72040209	\$ 384.91

#### Soft-Touch Sitters

Size	Lilac 	Aqua 	Dark Gray 	MSRP
Size 1	___ 77010102	___ 77010104	___ 77010109	\$ 740.64
Size 2	___ 77010202	___ 77010204	___ 77010209	\$ 974.55