


Please mark clearly with an X the items you wish to purchase.  
 Email form to Ki Mobility Canada Ltd. at [sales@kimobilitycanada.ca](mailto:sales@kimobilitycanada.ca) or fax form to 403-253-1744. Contact our customer service team at (877) 252-6655.

Account #: \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO #: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (     ) \_\_\_\_\_

### Special Tomato® MPS Features:

- Individualized adjustment allows for many years of growth
- Seamless latex-free cushions are impermeable to fluids and soft, yet durable
- Cleans easily with common disinfectants
- All Special Tomato MPS Seating Systems are FMVSS213 (Crash Testing) Compliant



### Select Your Cushions: Choose One of Each

MPS Seat Specifications	Small Seat Cushions				Large Seat Cushions			
	MS1	MS3	MS5	MS7	ML1	ML3	ML5	ML7
Target User Age	9 mos - 3 yrs	2 - 5 yrs	4 - 8 yrs	5 - 9 yrs	6 - 9 yrs	8 - 11 yrs	10 - 14 yrs	12 yrs - Adult
Target User Height	28" - 37"	32" - 43"	40" - 50"	43" - 52"	45" - 58"	48" - 60"	48" - 62"	58" - Adult
Maximum User Weight	80 lbs.				150 lbs.			
Seat Width	7"	9"	10"	10"	11"	12"	13"	13"
Seat Depth	6"	8"	10"	12"	11"	13"	15"	17"
Seat to Footrest	6" - 16"				8" - 18"			
Hip Flex Angle	-15° to +15°				-15° to +15°			
Floor to Seat Height	15"				18"			

MPS Back Specifications	Small Back Cushions		Large Back Cushions	
	MS2	MS4	ML2	ML4
Target User Age	9 mos - 3 yrs	2 - 9 yrs	6 - 9 yrs	8 yrs - Adult
Target User Height	28" - 37"	32" - 52"	45" - 58"	48" - Adult
Maximum User Weight	80 lbs.		150 lbs.	
Back Width	8"	9"	10"	11"
Back Height	8.5"	11"	11"	13"
Top of Back Cushion	11" - 16"	15" - 20"	15" - 20"	19" - 22"
Top of Head Cushion	18" - 23"	21" - 26"	24" - 28"	26" - 30"

Product	Small Dimensions	Large Dimensions	Small Weight	Large Weight
Tray	24"W × 17"D × 2"H		8 lbs.	
MPS Seat Only	15"W × 19"D × 42"H	19"W × 23"D × 47"H	21 lbs.	26 lbs.
MPS Floor Kit	19"W × 37"D × 37"H	22"W × 41"D × 42"H	39 lbs.	44 lbs.
MPS Push Chair Kit	25"W × 38"D × 42"H	25"W × 42"D × 47"H	68 lbs.	77 lbs.

\* ML2 designed to fit both the Small and Large system, providing added wider Back Cushion option for Small size systems.

Please mark clearly with an X the items you wish to purchase.  
Email form to Ki Mobility Canada Ltd. at [sales@kimobilitycanada.ca](mailto:sales@kimobilitycanada.ca) or fax form to 403-253-1744. Contact our customer service team at (877) 252-6655.

## Customize Your Size



**Small Shell (MS)**



**Large Shell (ML)**

**Shells**  
Choose a Shell size to determine your options

### Options for (MS) Shell

#### Headrests

MHS (Standard)



7.5"W x 6.5"H


MHL (w/ Laterals)



7"W x 6"H

#### Back Cushions

MS2



8"W x 8.5"H

MS4



9"W x 11"H

ML2



10"W x 11"H

#### Seat Cushions

MS1




7"W x 6"D

MS3



9"W x 8"D

MS5



10"W x 10"D

MS7



10"W x 12"D

### Options for (ML) Shell

#### Headrests

MHS (Standard)



7.5"W x 6.5"H

MHL (w/ Laterals)



7"W x 6"H


#### Back Cushions

ML2



10"W x 11"H


ML4



11"W x 13"H


#### Seat Cushions

ML1




11"W x 11"D

ML3



12"W x 13"D

ML5



13"W x 15"D

ML7



13"W x 17"D

Toll-Free: (877) 252-6655

Fax: (403) 253-1744

Email: [sales@kimobilitycanada.ca](mailto:sales@kimobilitycanada.ca)

Order Forms Available Online at: [www.kimobilitycanada.ca](http://www.kimobilitycanada.ca)

Please mark clearly with an X the items you wish to purchase.  
 Email form to Ki Mobility Canada Ltd. at sales@kimobilitycanada.ca or fax form to 403-253-1744. Contact our customer service team at (877) 252-6655.

## Customize Your Size



**MPS Shell with Footrest**

\_\_\_ Small - 79030109  
 \$ 2,623.16  
 \_\_\_ Large - 79030309  
 \$ 3,383.85



**MPS Mobile Floor Base**

\_\_\_ Small - 79010109  
 \$ 4,197.63  
 \_\_\_ Large - 79010309  
 \$ 5,422.58



**MPS Push Chair Base**

\_\_\_ Small - 79020109  
 \$ 6,669.12  
 \_\_\_ Large - 79020309  
 \$ 7,548.15

## Select MPS Cushions: One of Each

### Headrest Cushions

Size	Cherry	Lilac	Lime	Aqua	Dark Gray	Light Gray	MSRP
MHS	___ 79100051	___ 79100052	___ 79100053	___ 79100054	___ 79100059	___ 79100060	\$ 203.53
MHL	___ 79110051	___ 79110052	___ 79110053	___ 79110054	___ 79110059	___ 79110060	\$ 253.60

### Back Cushions

Size	Cherry	Lilac	Lime	Aqua	Dark Gray	Light Gray	MSRP
MS2	___ 79120551	___ 79120552	___ 79120553	___ 79120554	___ 79120559	___ 79120560	\$ 285.35
MS4	___ 79120151	___ 79120152	___ 79120153	___ 79120154	___ 79120159	___ 79120160	\$ 310.52
ML2	___ 79120751	___ 79120752	___ 79120753	___ 79120754	___ 79120759	___ 79120760	\$ 355.52
ML4	___ 79120351	___ 79120352	___ 79120353	___ 79120354	___ 79120359	___ 79120360	\$ 365.05

\* ML2 designed to fit both the Small and Large system, providing added wider Back Cushion option for Small size systems.

### Seat Cushions

Size	Cherry	Lilac	Lime	Aqua	Dark Gray	Light Gray	MSRP
MS1	___ 79050551	___ 79050552	___ 79050553	___ 79050554	___ 79050559	___ 79050560	\$ 285.34
MS3	___ 79050151	___ 79050152	___ 79050153	___ 79050154	___ 79050159	___ 79050160	\$ 310.52
MS5	___ 79050251	___ 79050252	___ 79050253	___ 79050254	___ 79050259	___ 79050260	\$ 335.70
MS7	___ 79050651	___ 79050652	___ 79050653	___ 79050654	___ 79050659	___ 79050660	\$ 355.51
ML1	___ 79050751	___ 79050752	___ 79050753	___ 79050754	___ 79050759	___ 79050760	\$ 355.51
ML3	___ 79050351	___ 79050352	___ 79050353	___ 79050354	___ 79050359	___ 79050360	\$ 365.05
ML5	___ 79050451	___ 79050452	___ 79050453	___ 79050454	___ 79050459	___ 79050460	\$ 374.59
ML7	___ 79050851	___ 79050852	___ 79050853	___ 79050854	___ 79050859	___ 79050860	\$ 396.69

\* MS1 Cushion is not certified for use in a motor vehicle.

\$

Toll-Free: (877) 252-6655      Fax: (403) 253-1744      Email: sales@kimobilitycanada.ca      **Your MSRP Total**  
 Order Forms Available Online at: www.kimobilitycanada.ca

Please mark clearly with an X the items you wish to purchase.  
 Email form to Ki Mobility Canada Ltd. at [sales@kimobilitycanada.ca](mailto:sales@kimobilitycanada.ca) or fax form to 403-253-1744. Contact our customer service team at (877) 252-6655.

## MPS Accessories



### Mobile Floor Base

Built-in tray storage on back. Four Swivel casters with rear locks. 10° - 25° of tilt.

Size	Part #	MSRP
Small	79070100	\$ 1,808.74
Large	79070300	\$ 2,315.10



### MPS Tray

Depth adjustable. 60° tilt adjustment. Tray easily attaches to back of Mobile Floor Base for storage.

Size	Part #	MSRP
Small	79060159	\$ 553.85
Large	79060359	\$ 624.60

### Push Chair Accessories

	Part #	MSRP
MPS Push Chair Base	79200000	\$ 5,040.28
MPS Sun Shade for Push Chair Base	79210000	\$ 831.66

### Push Chair w/ MPS Seat and Cushions



MPS Sun Shade for Push Chair Base