

Please mark clearly with an X the items you wish to purchase.
 Email form to NuVision Rehab Group at sales@nuvisionrehab.com or fax form to 403-253-1744. Contact our customer service team at (877) 252-6655.

Account #: _____ Ship to: _____
 Date: _____ Name: _____
 PO #: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

Soft-Touch® Sitters Features:

- Contoured Head & Lateral Support
- 5-point harness
- Attachment straps
- Anti-thrust seat surface
- Up to 25° of tilt



Mobile Base



Stationary Base



Floor Base

Sitter Only

Size	Aqua	Dark Gray	Lilac	MSRP
Size 1	77010104	77010109	77010102	\$ 660.10
Size 2	77010204	77010209	77010202	\$ 868.58
Size 3	77010304	77010309	77010302	\$ 1,104.86
Size 4	77010404	77010409	77010402	\$ 1,222.95
Size 5	77010504	77010509	77010502	\$ 1,862.17

Sitter on Floor Base

Size	Aqua	Dark Gray	Lilac	MSRP
Size 1	78010104	78010109	78010102	\$ 850.62
Size 2	78010204	78010209	78010202	\$ 1,005.33
Size 3	78010304	78010309	78010302	\$ 1,235.23

Optional Accessories

Item	Part #	MSRP
Floor Sitter Wedge (Fits Sizes 1-3)	78000009	\$ 542.29
Small Stationary Base (Fits Sizes 1-3)	77100100	\$ 953.12
Large Stationary Base (Fits Sizes 4-5)	77100200	\$ 1,174.77
Small Mobile Base (Fits Sizes 1-3)	77200100	\$ 1,168.64
Large Mobile Base (Fits Sizes 4-5)	77200200	\$ 1,429.38
Headrest	77020000	\$ 253.32

Sitter on Stationary Base

Size	Aqua	Dark Gray	Lilac	MSRP
Size 1	77110104	77110109	77110102	\$ 1,185.48
Size 2	77110204	77110209	77110202	\$ 1,420.93
Size 3	77110304	77110309	77110302	\$ 1,715.17
Size 4	77110404	77110409	77110402	\$ 1,727.55
Size 5	77110504	77110509	77110502	\$ 2,705.89

Sitter on Mobile Base

Size	Aqua	Dark Gray	Lilac	MSRP
Size 1	77210104	77210109	77210102	\$ 1,612.83
Size 2	77210204	77210209	77210202	\$ 1,796.39
Size 3	77210304	77210309	77210302	\$ 2,086.77
Size 4	77210404	77210409	77210402	\$ 2,300.68
Size 5	77210504	77210509	77210502	\$ 3,240.77



Soft-Touch® Sitter on Chair



Optional Headrest



5-Point Harness

Please mark clearly with an X the items you wish to purchase.
 Email form to NuVision Rehab Group at sales@nuvisionrehab.com or fax form to 403-253-1744. Contact our customer service team at (877) 252-6655.

Base Sizing Chart

Sitter Base Specifications	Floor Base	Stationary Base		Mobile Base	
		Small	Large	Small	Large
Use With Sitter Sizes	1, 2, 3	1, 2, 3	4, 5	1, 2, 3	4, 5
Overall Width	19"	17.5"	19.25"	17.5"	19.25"
Overall Depth	22.5"	30"	35"	30"	35"
Overall Height	12.5"	16"	23"	19.5"	26.5"
Weight of Base	14 lbs.	12 lbs.	15 lbs.	14 lbs.	17 lbs.
Degree of Tilt	25°	25°		25°	

Sitter Sizing Chart

Dimensions	Size 1	Size 2	Size 3	Size 4	Size 5
Target User Weight	20 - 30 lbs.	30 - 60 lbs.	40 - 80 lbs.	50 - 110 lbs.	90 - 200 lbs.
Target User Height	30" - 35"	35" - 48"	40" - 56"	50" - 65"	54" - 74"
Target User Age	1 - 2 yrs	2 - 6 yrs	6 - 9 yrs	9 - 14 yrs	13 yrs - Adult
Inside Seat Depth	7"	9.5"	11.5"	13.5"	16"
Inside Seat Width	8"	10"	12"	14"	16"
Inside Seat Height	21"	26"	29"	35"	39"
Floor To Seat Ht. Upright on Stationary Base	10"	10"	9.75"	11.5"	12.5"
Floor To Seat Ht. Tilted on Stationary Base	10.5"	10.5"	10.5"	13"	14.25"
Floor To Seat Ht. Upright on Mobile Base	12.5"	13.25"	13"	14.75"	16"
Floor To Seat Ht. Tilted on Mobile Base	12.75"	13.75"	13.75"	16"	17.75"
Floor To Seat Ht. Upright on Floor Base	5"	5.5"	5.5"	DNA	DNA
Floor To Seat Ht. Tilted on Floor Base	4.5"	4.5"	5"	DNA	DNA
Overall Seat Depth	12"	14"	16"	19"	21"
Overall Seat Depth	12"	14"	16"	18"	20"
Overall Seat Depth	25"	29"	34"	40"	44"
Weight of Sitter	6 lbs.	10 lbs.	13 lbs.	15 lbs.	19 lbs.