

Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Effective 05/10/2018

Liberty FT **\$3725**
 250 lb weight capacity. Tilt range: 0° - 20°.

Frame

Seat Width Range / Seat Width

Select Seat Width within desired Seat Width Range

Seat Width Range	Seat Width				
Small 14" - 18"	<input type="checkbox"/> 14 N/C	<input type="checkbox"/> 15 N/C	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C
Medium 16" - 20"	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 N/C
Large 18" - 22"	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 N/C	<input type="checkbox"/> 21 N/C	<input type="checkbox"/> 22 N/C

Seat Depth Range / Seat Depth

Select Seat Depth within desired Seat Depth Range

Seat Depth Range	Seat Depth		
Short 14" - 16"	<input type="checkbox"/> 14 N/C	<input type="checkbox"/> 15 N/C	<input type="checkbox"/> 16 N/C
Medium 16" - 18"	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C
Long 18" - 20"	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 N/C

Assistive Gas Strut Set-Up		
Client Weight	Setting	Spring Type (Newtons)
<input type="checkbox"/> > 200 lbs	1	450 N
<input type="checkbox"/> 176 - 200 lbs	2	450 N
<input type="checkbox"/> 151 - 175 lbs ¹	3	450 N
<input type="checkbox"/> 151 - 175 lbs ²	1	300 N
<input type="checkbox"/> 125 - 150 lbs	2	300 N
<input type="checkbox"/> < 125 lbs	3	300 N

¹ Standard on chairs less than 16" deep
² Standard on chairs 16" deep or larger

Seat Option

- Standard Seat Sling _____ Std
- Omit Seat _____ N/C

Transit Option

- Transit System _____ \$320
Available only with Fold Down Back

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range

Back Post	Price	Height Range	Back Height				
Stroller Handle Back Post ¹	Std	Short	<input type="checkbox"/> 20				
		Tall	<input type="checkbox"/> 24				
Height Adjustable Straight with Push Handle	\$36	Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20
		Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22
Height Adjustable 8° Bend with Push Handle ²	\$36	Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20
		Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22

¹ Not available with Back Upholstery

² Not available with IV Pole

Stroller Handle Back Post Option

Available only with Stroller Handle Back Post

Removable Stroller Handle Extension _____ \$140

Backrest Options

Fold Down Back _____ \$95

Back Upholstery

Standard _____ Std

Tension Adjustable _____ \$295
 Not available on Seat Widths less than 16" combined with 21" and 22" Backrest Heights

Omit Back Upholstery _____ - \$20

Footrest

Hanger

60° _____ Std
 Available only with Extension Tube Hangers

70° _____ N/C

80° _____ N/C

90° _____ N/C

Available only with Front Mount Hangers. Not available with 8" Casters combined with 3" Extension

Elevating Leg Rests (ELR) _____ \$369

Pro Elevating Leg Rests (PRO ELR) _____ \$500

Omit Hangers _____ - \$25

Footplate

Composite _____ Std
 Not available with Front Mount Hangers or PRO ELR

Aluminum Angle Adjustable _____ \$165

One Piece Flip Up Angle Adjustable _____ \$238
 Not available with Front Mount Hangers, ELR or PRO ELR

Footrest continued

Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length												
Extension Tube Hangers (60°, 70°, 80°)	Short	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5					
	Medium	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5				
	Long	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18				
Front Mount Hangers ¹	Standard Selections										Requires 3" Extension			
	Short (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5						<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	
		<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8									<input type="checkbox"/> 10	<input type="checkbox"/> 10.5
	Medium (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5					<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5
		<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10						<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
	Long (70°, 80°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5		
<input type="checkbox"/> 9		<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12			<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15		
Elevating Leg Rests (ELR)	Short	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5								
	Long	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5							
Extended Length ELR	Short	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20.5								
	Long	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21.5	<input type="checkbox"/> 22.5						
Pro ELR Low Pivot Setting	Short	<input type="checkbox"/> 10.5 - 15.5												
	Long	<input type="checkbox"/> 13 - 19												
Pro ELR High Pivot Setting	Short	<input type="checkbox"/> 8.5 - 17												
	Long	<input type="checkbox"/> 11.5 - 17.5												

¹ See Minimum Seat Height chart below

Minimum Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

Heel Loop

- Heel Loops _____ \$45
Not available with One Piece Flip Up Angle Adjustable Footplate
- Omit Heel Loops _____ N/C

Calf Strap

- Velcro Adjustable _____ \$44
- Padded Velcro Adjustable _____ \$85
- Bodypoint Padded AeroMesh _____ \$123

3" Extension

Available only with Front Mount Hangers. Provided as parts.

- 3" Extension for Front Mount Hanger _____ N/C

Residual Limb Support

- Residual Limb Support
 - Left _____ \$450
 - Right _____ \$450
 - Omit Left Hanger _ -\$12
 - Omit Right Hanger _ -\$12

Cover

- Smooth
- Air Mesh


Pad Size

- | | | | |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 | <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

Footrest continued

Hanger Release

- 4-Way - In line _____ Std
 See reference A to the right. Not available with ELR or PRO ELR.
- 4-Way - Out _____ N/C
 See reference B to the right



Seat Height, Casters and Rear Wheels

Seat Height / Caster Size / Rear Wheel Size

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size			
	4"	5"	6"	8"
12"	Seat Height			
	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 17
	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 18
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 19
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	
	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	
		<input type="checkbox"/> 18	<input type="checkbox"/> 19	
16"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 17
	<input type="checkbox"/> 14 ¹	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 18
	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 19
	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 19.5
		<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 20
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	
		<input type="checkbox"/> 18	<input type="checkbox"/> 19	
		<input type="checkbox"/> 18.5		
20"	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 17
	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 18
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 19
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 19.5
	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 20
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	
		<input type="checkbox"/> 18	<input type="checkbox"/> 19	
		<input type="checkbox"/> 18.5		

Rear Wheel Size	Caster Size			
	4"	5"	6"	8"
22"	Seat Height			
	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 17
	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 18
	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 19
	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 19.5
		<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 20
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	
		<input type="checkbox"/> 18	<input type="checkbox"/> 19	
24"	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 17
	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 18
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 19
		<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 19.5
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 20
		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	
		<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	

¹ Available only with Low Profile tires

Seat Height, Casters and Rear Wheels continued

Caster Type

- 1" Poly _____ Std
- 1" Poly Aluminum _____ \$138
Not available with 6" and 8" Caster Sizes
- 1" Pneumatic _____ \$200
Not available with 4" and 5" Caster Sizes
- 1.5" Poly _____ \$90
Not available with 8" Caster Size
- 1.5" Soft Roll Aluminum _____ \$175
Not available with 8" Caster Size
- 2" Poly _____ \$63
Available only with 6" Caster Size
- 2" Pneumatic w/ Foam Insert _____ \$315
Available only with 8" Caster Size

Rear Wheel Type

- 24" Spoke _____ N/C
- Spoke (20" or 22") _____ \$180
- Mag _____ Std

Rear Tire Type

Full Profile

- Street Pneumatic _____ Std
- Full Poly _____ N/C
- Pneumatic with Airless Insert _____ \$80
- Euro Trek Knobby _____ \$119
Available only on 24" rear wheels

Low Profile

- Low Poly _____ N/C
Available only with 16" rear wheel
- Sentinel Solid Tire _____ \$120
Not available on 12" wheels, 16" wheels or 20" and 22" Mag wheels
- Schwalbe Marathon Plus Evolution _____ \$125
Available only on 24" spoke wheels

Wheel Locks

Wheel Lock

- Push to Lock _____ Std
- Pull to Lock _____ N/C
Not available with 12" or 16" Full Poly tires
- Omit Wheel Locks _____ N/C

Handrim

Handrims not available on 12" or 16" rear wheels

Standard Handrims

- Aluminum Anodized _____ Std
- Plastic Coated _____ \$95
- Projection (8 Vertical Projections) _____ \$299
Available on 22" and 24" Spoke wheels & 24" Mag wheels w/full tires
- Omit Handrims _____ N/C

Natural Fit Handrims

Not available on 20" rear wheels or 22" Mag wheels. Not recommended with Mag wheels.

- Natural Fit LT Regular Grip _____ \$445
- Natural Fit Standard Regular Grip _____ \$445
- Natural Fit LT Super Grip _____ \$495
- Natural Fit Standard Super Grip _____ \$495
- Natural Fit LT Omit Grip _____ \$345
- Natural Fit Standard Omit Grip _____ \$345

Handrim Option

Not available on Plastic Coated or Projection Handrims

- 3M Micro Grip _____ \$81

Axle

- Fixed _____ N/C
- Quick Release _____ Std

Wheel Lock Options

- Extension Handles _____ \$84
Available only with Push to Lock or Pull to Lock Wheel Locks
- Grade Aids _____ \$187
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires

Armrests

Armrest Please select Arm Pad option

- Height Adjustable T-Arm (8.5" - 13.5") \$285
- Angle Adjustable Locking Extendable Flip Up \$295
Available only with Stroller Handle Back Post
- Omit Armrests N/C

Arm Pad

- Standard Arm Pad N/C Desk Full
- Waterfall Arm Pad \$35 Desk Full
Not available with Angle Adjustable Locking Extendable Flip Up Armrests
- Foam Grip N/C Desk Full
Available only with Angle Adjustable Locking Extendable Flip Up Armrests

Armrest Continued

- Tubular Flip Up Std
Available only with Stroller Handle Back Post. Not available with Arm Pads.

Accessories

Removable Side Guard

- Composite Side Guards \$175
Not available with T-Arms

Anti-Tip

- Rear Anti-Tips \$125

Positioning Belt

- 1 1/2" Auto Buckle \$80
- 1 1/2" Padded Auto Buckle \$170
- 2" Airline Buckle \$150
- Bodypoint Evoflex - Pelvic Stabilizer \$219
Size: Small (16-24") Medium (19-29") Large (25-37")

Belt Mounting Kit

- Seat Belt Mounting Hardware \$38
- Bodypoint 1" Band Clamp Belt Mounting Kit \$126

O₂ Holder

- O₂ Holder \$365
Mounted to the left side. Not available with Fold Down Back.

Other Accessories

- Tool Kit \$60
- Removable Underseat Pouch \$75
- Ki Mobility Backpack \$90
- Black Neoprene Impact Guards \$95
- Clear Spoke Guards \$144
Not available with Mag wheels
- Black Spoke Guards \$144
Not available with 20" wheels or Mag wheels
- Cane and Crutch Holder \$195
- IV Pole \$175
Not available with Transit Option or 8° Bend with Push Handle Back Post. Left side mount not available with O₂ Holder.
 Right Side Left Side

Frame Color

Standard Colors

Gloss Finish

- Black N/C
- Candy Blue N/C
- Candy Purple N/C
- Candy Red N/C
- Copper N/C
- Silver N/C

Matte Finish

- Army Green N/C
- Charcoal Gray N/C
- Iridescent Burgundy N/C
- Iridescent Navy N/C
- Shadow Black N/C

Optional Colors

Extended lead time

Gloss Finish

- Blue Skies \$90
- Bubblegum Pink \$90
- Fuchsia \$90
- Granny Smith \$90
- Metallic Orange \$90
- Teal \$90
- Yellow \$90

Matte Finish

- Emerald Green \$90
- Retro Red \$90
- Royal Blue \$90
- Smooth White \$90

AXIOM

Axiom G, S and P Cushions

Select size in chart at right. Weight capacity: 350 lbs.

Cushions Under 22" Wide

- Axiom G _____ \$195
- Axiom S _____ \$270
- Axiom P _____ \$575

22" Wide Cushions

- Axiom G _____ \$275
- Axiom S _____ \$400
- Axiom P _____ \$695

Axiom G, S or P Cushion Quantity: 1 2

If no selection is made, quantity 1 will be selected

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Depth	14"	<input type="checkbox"/>								
	15"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	17"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	18"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>
	19"					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	20"			<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>

Axiom SP Cushions

Select size in chart at right. Weight capacity: 350 lbs.

Cushions Under 22" Wide

- Axiom SP - Visco _____ \$640
- Axiom SP - Fluid _____ \$695

22" Wide Cushions

- Axiom SP - Visco _____ \$775
- Axiom SP - Fluid _____ \$830

Axiom SP Cushion Quantity: 1 2

If no selection is made, quantity 1 will be selected

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Depth	14"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
	15"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>
	17"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>
	18"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Axiom Accessories

Outer Covers

- Axiom G Extra Outer Cover _____ \$95
- Axiom S Extra Outer Cover _____ \$95
- Axiom P Extra Outer Cover _____ \$95
- Axiom SP Extra Visco Outer Cover _____ \$134
- Axiom SP Extra Fluid Outer Cover _____ \$134

Outer Cover Quantity: 1 2

If no selection is made, quantity 1 will be selected

Inner Covers

- Axiom SP Extra Visco Inner Cover _____ \$134
- Axiom SP Extra Fluid Inner Cover _____ \$134

Inner Cover Quantity: 1 2

If no selection is made, quantity 1 will be selected

Solid Seat Insert

Select size in chart below. Weight capacity: 350 lbs.

Solid Seat Insert _____ \$73

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Depth	14"	<input type="checkbox"/>								
	15"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	17"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	18"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>
	19"					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	20"			<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>

Solid Seat Insert Quantity: 1 2

If no selection is made, quantity 1 will be selected



Elite Back

Weight Capacity: 300 lbs.

STEP 1: Choose Back Size

Quick Ship

Quick Ship products are stocked

- Matrx Elite Back 14 x 14 - PBE1414 \$766
- Matrx Elite Back 14 x 16 - PBE1416 \$766
- Matrx Elite Back 14 x 18 - PBE1418 \$766
- Matrx Elite Back 15 x 14 - PBE1514 \$766
- Matrx Elite Back 15 x 16 - PBE1516 \$766
- Matrx Elite Back 15 x 18 - PBE1518 \$766
- Matrx Elite Back 16 x 12 - PBE1612 \$766
- Matrx Elite Back 16 x 14 - PBE1614 \$766
- Matrx Elite Back 16 x 16 - PBE1616 \$766
- Matrx Elite Back 16 x 18 - PBE1618 \$766
- Matrx Elite Back 16 x 20 - PBE1620 \$766
- Matrx Elite Back 17 x 16 - PBE1716 \$766
- Matrx Elite Back 17 x 18 - PBE1718 \$766
- Matrx Elite Back 17 x 20 - PBE1720 \$766
- Matrx Elite Back 18 x 14 - PBE1814 \$766
- Matrx Elite Back 18 x 16 - PBE1816 \$766
- Matrx Elite Back 18 x 18 - PBE1818 \$766
- Matrx Elite Back 18 x 20 - PBE1820 \$766
- Matrx Elite Back 19 x 16 - PBE1916 \$766
- Matrx Elite Back 19 x 18 - PBE1918 \$766
- Matrx Elite Back 20 x 16 - PBE2016 \$766
- Matrx Elite Back 20 x 18 - PBE2018 \$766
- Matrx Elite Back 20 x 20 - PBE2020 \$766

On Demand

On Demand products are ordered with an extended lead time

- Matrx Elite Back 12 x 10 - PBE1210 \$766
- Matrx Elite Back 12 x 12 - PBE1212 \$766
- Matrx Elite Back 12 x 14 - PBE1214 \$766
- Matrx Elite Back 12 x 16 - PBE1216 \$766
- Matrx Elite Back 12 x 18 - PBE1218 \$766
- Matrx Elite Back 13 x 10 - PBE1310 \$766
- Matrx Elite Back 13 x 12 - PBE1312 \$766
- Matrx Elite Back 13 x 14 - PBE1314 \$766
- Matrx Elite Back 13 x 16 - PBE1316 \$766
- Matrx Elite Back 13 x 18 - PBE1318 \$766
- Matrx Elite Back 14 x 10 - PBE1410 \$766
- Matrx Elite Back 14 x 12 - PBE1412 \$766
- Matrx Elite Back 14 x 20 - PBE1420 \$766
- Matrx Elite Back 15 x 10 - PBE1510 \$766
- Matrx Elite Back 15 x 12 - PBE1512 \$766
- Matrx Elite Back 15 x 20 - PBE1520 \$766
- Matrx Elite Back 16 x 10 - PBE1610 \$766
- Matrx Elite Back 17 x 10 - PBE1710 \$766
- Matrx Elite Back 17 x 12 - PBE1712 \$766
- Matrx Elite Back 17 x 14 - PBE1714 \$766
- Matrx Elite Back 18 x 10 - PBE1810 \$766
- Matrx Elite Back 18 x 12 - PBE1812 \$766
- Matrx Elite Back 19 x 10 - PBE1910 \$766
- Matrx Elite Back 19 x 12 - PBE1912 \$766
- Matrx Elite Back 19 x 20 - PBE1920 \$766
- Matrx Elite Back 20 x 10 - PBE2010 \$766
- Matrx Elite Back 20 x 12 - PBE2012 \$766

STEP 2: Choose Back Hardware (Selection Required)

Back Hardware

- Easy Set Hardware - PBMH _____ N/C
*Not available with Angle Adjustable Locking Flip Up or Tubular Flip Up Armrests,
 O2 Holder or welded rigidizer bar*
- Split Mounting Hardware - SMHU _____ N/C



Elite Deep Back

Weight Capacity: 300 lbs.

STEP 1: Choose Back Size

Quick Ship

Quick Ship products are stocked

<input type="checkbox"/>	Matrx Elite Deep Back 14 x 16 - EDB1416	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 14 x 18 - EDB1418	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 15 x 16 - EDB1516	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 15 x 18 - EDB1518	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 16 x 16 - EDB1616	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 16 x 18 - EDB1618	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 16 x 20 - EDB1620	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 17 x 16 - EDB1716	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 17 x 18 - EDB1718	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 17 x 20 - EDB1720	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 18 x 16 - EDB1816	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 18 x 18 - EDB1818	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 18 x 20 - EDB1820	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 19 x 18 - EDB1918	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 19 x 20 - EDB1920	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 20 x 16 - EDB2016	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 20 x 18 - EDB2018	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 20 x 20 - EDB2020	\$918

On Demand

On Demand products are ordered with an extended lead time

<input type="checkbox"/>	Matrx Elite Deep Back 14 x 12 - EDB1412	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 14 x 20 - EDB1420	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 15 x 10 - EDB1510	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 15 x 12 - EDB1512	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 15 x 20 - EDB1520	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 16 x 10 - EDB1610	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 16 x 12 - EDB1612	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 17 x 10 - EDB1710	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 17 x 12 - EDB1712	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 18 x 10 - EDB1810	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 18 x 12 - EDB1812	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 19 x 10 - EDB1910	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 19 x 12 - EDB1912	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 19 x 16 - EDB1916	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 20 x 10 - EDB2010	\$918
<input type="checkbox"/>	Matrx Elite Deep Back 20 x 12 - EDB2012	\$918

STEP 2: Choose Back Hardware (Selection Required)

Back Hardware

- Easy Set Hardware - PBMH _____ N/C
Not available with Angle Adjustable Locking Flip Up or Tubular Flip Up Armrests,
 O2 Holder or welded rigidizer bar
- Split Mounting Hardware - SMHU _____ N/C



Elite TR Back

Weight Capacity: 300 lbs.

STEP 1: Choose Back Size

Quick Ship

Quick Ship products are stocked

- Matrx Elite Tilt / Recline Back 15 x 18 - PBE1518-TR ____ \$866
- Matrx Elite Tilt / Recline Back 16 x 18 - PBE1618-TR ____ \$866
- Matrx Elite Tilt / Recline Back 16 x 20 - PBE1620-TR ____ \$866
- Matrx Elite Tilt / Recline Back 17 x 18 - PBE1718-TR ____ \$866
- Matrx Elite Tilt / Recline Back 17 x 20 - PBE1720-TR ____ \$866
- Matrx Elite Tilt / Recline Back 18 x 18 - PBE1818-TR ____ \$866
- Matrx Elite Tilt / Recline Back 18 x 20 - PBE1820-TR ____ \$866
- Matrx Elite Tilt / Recline Back 19 x 18 - PBE1918-TR ____ \$866
- Matrx Elite Tilt / Recline Back 19 x 20 - PBE1920-TR ____ \$866
- Matrx Elite Tilt / Recline Back 20 x 18 - PBE2018-TR ____ \$866
- Matrx Elite Tilt / Recline Back 20 x 20 - PBE2020-TR ____ \$866

On Demand

On Demand products are ordered with an extended lead time

- Matrx Elite Tilt / Recline Back 15 x 16 - PBE1516-TR ____ \$866
- Matrx Elite Tilt / Recline Back 15 x 20 - PBE1520-TR ____ \$866
- Matrx Elite Tilt / Recline Back 16 x 16 - PBE1616-TR ____ \$866
- Matrx Elite Tilt / Recline Back 17 x 16 - PBE1716-TR ____ \$866
- Matrx Elite Tilt / Recline Back 18 x 16 - PBE1816-TR ____ \$866
- Matrx Elite Tilt / Recline Back 19 x 16 - PBE1916-TR ____ \$866
- Matrx Elite Tilt / Recline Back 20 x 16 - PBE2016-TR ____ \$866

STEP 2: Choose Back Hardware (Selection Required)

Back Hardware

- Easy Set Hardware - PBMH _____ N/C
*Not available with Angle Adjustable Locking Flip Up or Tubular Flip Up Armrests,
 O2 Holder or welded rigidizer bar*
- Split Mounting Hardware - SMHU _____ N/C



Matrix Elan Headrest

Quick Ship

Quick Ship products are stocked

- Elan Hardware - MEHW _____ \$266
- Elan Mini Hardware - MEMH _____ \$254
- Elan Standard Pad 6" W x 3.5" H w/ Std Cover - ESP6 _____ \$310
- Elan Standard Pad 10" W x 5" H w/ Std Cover - ESP10 _____ \$310
- Elan Standard Pad 14" W x 5" H w/ Std Cover - ESP14 _____ \$310

On Demand

On Demand products are ordered with an extended lead time

- Elan Standard Pad 6" W x 3.5" H w/ Infection Control Cover - ESP6-IC _____ \$335
- Elan Standard Pad 10" W x 5" H w/ Infection Control Cover - ESP10-IC _____ \$335
- Elan Standard Pad 14" W x 5" H w/ Infection Control Cover - ESP14-IC _____ \$335
- Elan Occipital Pad 9" W x 7" H w/ Std Cover - EOP9 _____ \$341
- Elan Occipital Pad 9" W x 7" H w/ Infection Control Cover - EOP9-IC _____ \$366
- Elan Occipital Pad 12" W x 8" H w/ Std Cover - EOP12 _____ \$341
- Elan Occipital Pad 12" W x 8" H w/ Infection Control Cover - EOP12-IC _____ \$366
- Elan 4 Point Pad 11" W x 10" H w/ Std Cover - E4POINT _____ \$448
- Elan 4 Point Pad 11" W x 10" H w/ Infection Control Cover - E4POINT-IC _____ \$473

Matrix Elite/Mini Laterals

Quick Ship

Quick Ship products are stocked

- Elite Swing Away Lateral Hardware Only - SAHEB
 - Quantity 1 _____ \$300
 - Quantity 2 _____ \$600
- Elite Offset Swing Away Lateral Hardware Only (Left) - OSAH-L _____ \$300
- Elite Offset Swing Away Lateral Hardware Only (Right) - OSAH-R _____ \$300
- Pad for Swing Away Lateral Medium 3.75" H x 5.75" L - EPAD-M
 - Quantity 1 _____ \$147
 - Quantity 2 _____ \$294
- Pad for Swing Away Lateral Large 4.25" H x 6.5" L - EPAD-L
 - Quantity 1 _____ \$147
 - Quantity 2 _____ \$294

On Demand

On Demand products are ordered with an extended lead time

- Elite Lateral Fixed Standard Pad 5" H x 7.5" L (Left) - PBLS-L _____ \$116
- Elite Lateral Fixed Standard Pad 5" H x 7.5" L (Right) - PBLS-R _____ \$116
- Elite Lateral Fixed Small Pad 3" H x 5.5" L (Left) - PBSS-L _____ \$113
- Elite Lateral Fixed Small Pad 3" H x 5.5" L (Right) - PBSS-R _____ \$113
- Elite Lateral Fixed Offset Pad 5" H x 7.5" L (Left) - PBOL-L _____ \$125
- Elite Lateral Fixed Offset Pad 5" H x 7.5" L (Right) - PBOL-R _____ \$125
- Mini Swing-Away Lateral Hardware Only - SAHMINI
 - Quantity 1 _____ \$300
 - Quantity 2 _____ \$600
- Pad for Swing Away Lateral X Small 3.25" H x 4.25" L - EPAD-XS
 - Quantity 1 _____ \$147
 - Quantity 2 _____ \$294
- Pad for Swing Away Lateral Small 3.5" H x 5.25" L - EPAD-S
 - Quantity 1 _____ \$147
 - Quantity 2 _____ \$294

Hardware Options

On Demand

On Demand products are ordered with an extended lead time

- Easy Set Hardware (Pair) - PBMH _____ \$291
- Split Mounting Hardware (Pair) - SMH _____ \$291
- Heavy Duty Easy Set Hardware (Pair) - HDMH _____ \$413
- Extended Range Mounting Hardware (Pair) - ERMH _____ \$225
- Fixed Space Saver Mounting Hardware (Pair) - FSSH _____ \$233
- Carbon Fiber Fixed Mount Hardware for MX2 Back - IMX2CFH _____ \$325
- Quick Release Hardware for MX2 Back - IMX2QRH _____ \$360

Other Matrix Products

On Demand

On Demand products are ordered with an extended lead time

- Chest Harness Interface Kit for Elite Back - CHIEB _____ \$103
- Chest Harness Interface Kit for Mini Back - HIMINI _____ \$114
- Extended Mounting Pins for Mounting to Backs 2" Wider - EMP _____ \$60
- Privacy Flap X Small (Fits 12" to 14" Backs) - PFLAP-XS _____ \$53
- Privacy Flap Small (Fits 15" to 16" Backs) - PFLAP-S _____ \$53
- Privacy Flap Medium (Fits 17" to 18" Backs) - PFLAP-M _____ \$53
- Privacy Flap Large (Fits 19" to 20" Backs) - PFLAP-L _____ \$53
- Privacy Flap X Large (Fits 21" to 22" Backs) - PFLAP-XL _____ \$53