

Please mark clearly with an X the items you wish to purchase.  
 Email form to NuVision Rehab Group at [sales@nuvisionrehab.com](mailto:sales@nuvisionrehab.com) or fax form to 403-253-1744.  
 Contact our customer service team at (877) 252-6655.

Account #: \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO #: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (     ) \_\_\_\_\_

### Special Tomato Hi-Low™ Features:

- Affordable Customized Seating
- Tilt-in-Space with Child in Seat
- Floor to Table Adjustment
- Contoured Cushion Support
- Tool Free Adjustments
- Ergo Harness Design
- Soft-Touch® Material



Shown here with Standard Headrest (sold separately)

Adjustable 5-point harness

Shown here with Tray (sold separately)

Rear-locking casters

Adjustable Trunk cushions

Hip Flex™ Seat Cushion adjusts from 15° of hip flexion to 15° of hip extension



## Select Your Cushions: Choose One of Each

MPS Seat Specifications	Small Seat Cushions				Large Seat Cushions			
	MS1	MS3	MS5	MS7	ML1	ML3	ML5	ML7
Target User Age	9 mos - 3 yrs	2 - 5 yrs	4 - 8 yrs	5 - 9 yrs	6 - 9 yrs	8 - 11 yrs	10 - 14 yrs	12 yrs - Adult
Target User Height	28" - 37"	32" - 43"	40" - 50"	43" - 52"	45" - 58"	48" - 60"	48" - 62"	58" - Adult
Maximum User Weight	80 lbs.				130 lbs.			
Seat Width	7"	9"	10"	10"	11"	12"	13"	13"
Seat Depth	6"	8"	10"	12"	11"	13"	15"	17"
Seat to Footrest	6" - 16"				8" - 18"			
Hip Flex Angle	-15° to +15°				-15° to +15°			
Seat Tilt Angle Forward	10°				10°			
Seat Tilt Angle Back	-30°				-30°			
Floor to Seat Height	14" - 28.5"				17" - 28.5"			

MPS Back Specifications	Small Back Cushions		Large Back Cushions	
	MS2	MS4	ML2	ML4
Target User Age	9 mos - 3 yrs	2 - 9 yrs	6 - 9 yrs	8 yrs - Adult
Target User Height	28" - 37"	32" - 52"	45" - 58"	48" - Adult
Maximum User Weight	80 lbs.		130 lbs.	
Back Width	8"	9"	10"	11"
Back Height	8.5"	11"	11"	13"
Top of Back Cushion	11" - 16"	15" - 20"	15" - 20"	19" - 22"
Top of Head Cushion	18" - 23"	21" - 26"	24" - 28"	26" - 30"

Tray Specifications	
Dimensions	24"W x 17"D x 2"H
Weight	8 lbs.

Please mark clearly with an X the items you wish to purchase.  
Email form to NuVision Rehab Group at [sales@nuvisionrehab.com](mailto:sales@nuvisionrehab.com) or fax form to 403-253-1744.  
Contact our customer service team at (877) 252-6655.

## Customize Your Size



Small Shell (MS)



Large Shell (ML)

**Shells**  
Choose a Shell size to determine your options

### Options for (MS) Shell

#### Headrests

MHS (Standard)



7.5"W x 6.5"H

MHL (w/ Laterals)



7"W x 6"H

#### Back Cushions

MS2



8"W x 8.5"H

MS4



9"W x 11"H

ML2



10"W x 11"H


#### Seat Cushions

MS3



9"W x 8"D

MS5



10"W x 10"D

MS7



10"W x 12"D

### Options for (ML) Shell

#### Headrests

MHS (Standard)



7.5"W x 6.5"H

MHL (w/ Laterals)



7"W x 6"H


#### Back Cushions

ML2



10"W x 11"H


ML4



11"W x 13"H


#### Seat Cushions

ML1




11"W x 11"D

ML3




12"W x 13"D

ML5



13"W x 15"D

ML7



13"W x 17"D

Please mark clearly with an X the items you wish to purchase.  
 Email form to NuVision Rehab Group at sales@nuvisionrehab.com or fax form to 403-253-1744.  
 Contact our customer service team at (877) 252-6655.



### MPS Hi-Low™ Seating System

Includes Seat Shell with Footrest, 5-Point Harness, Hi-Low™ Base with Casters (Cushions Separate)

Size	Part #	MSRP
Small	79300109	\$ 4,511.51
Large	79300309	\$ 5,868.32



### Optional MPS Tray

Sold Separately

Size	Part #	MSRP
Small	79060159	\$ 553.85
Large	79060359	\$ 624.60

## Select MPS Cushions: One of Each



### Headrest Cushions

Size	Cherry	Lilac	Lime	Aqua	Dark Gray	Light Gray	MSRP
MHS	79100051	79100052	79100053	79100054	79100059	79100060	\$ 203.53
MHL	79110051	79110052	79110053	79110054	79110059	79110060	\$ 253.60

### Back Cushions

Size	Cherry	Lilac	Lime	Aqua	Dark Gray	Light Gray	MSRP
MS2	79120551	79120552	79120553	79120554	79120559	79120560	\$ 285.35
MS4	79120151	79120152	79120153	79120154	79120159	79120160	\$ 310.52
ML2	79120751	79120752	79120753	79120754	79120759	79120760	\$ 355.52
ML4	79120351	79120352	79120353	79120354	79120359	79120360	\$ 365.05

### Seat Cushions

Size	Cherry	Lilac	Lime	Aqua	Dark Gray	Light Gray	MSRP
MS1	79050551	79050552	79050553	79050554	79050559	79050560	\$ 285.34
MS3	79050151	79050152	79050153	79050154	79050159	79050160	\$ 310.52
MS5	79050251	79050252	79050253	79050254	79050259	79050260	\$ 335.70
MS7	79050651	79050652	79050653	79050654	79050659	79050660	\$ 355.70
ML1	79050751	79050752	79050753	79050754	79050759	79050760	\$ 355.70
ML3	79050351	79050352	79050353	79050354	79050359	79050360	\$ 365.05
ML5	79050451	79050452	79050453	79050454	79050459	79050460	\$ 374.59
ML7	79050851	79050852	79050853	79050854	79050859	79050860	\$ 396.69

Does the intended user of this Hi-Low system weigh over 95 lbs. (43kg)?  Yes  No

\$

### \* Warranty:

5-Year Manufacturer's Warranty on cushions and seat shell  
 2-Year Manufacturer's Warranty on Hi-Low™ base and strapping

**Your MSRP Total**