



Soft-Touch® Sitters

30DPDFODDLDLPRXLBSXFD
 PDLOIRPRLRLLOLDGDDDO#LPRLOLFDDGDFDRDIRPR
 RDRFXFRPFFDPD

Base Sizing Chart

Sitter Base Specifications	Floor Base	Stationary Base		Mobile Base	
		Small	Large	Small	Large
Use With Sitter Sizes	1, 2, 3	1, 2, 3	4, 5	1, 2, 3	4, 5
Overall Width	19"	17.5"	19.25"	17.5"	19.25"
Overall Depth	22.5"	30"	35"	30"	35"
Overall Height	12.5"	16"	23"	19.5"	26.5"
Weight of Base	14 lbs.	12 lbs.	15 lbs.	14 lbs.	17 lbs.
Degree of Tilt	25°	25°		25°	

Sitter Sizing Chart

Dimensions	Size 1	Size 2	Size 3	Size 4	Size 5
Target User Weight	20 - 30 lbs.	30 - 60 lbs.	40 - 80 lbs.	50 - 110 lbs.	90 - 200 lbs.
Target User Height	30" - 35"	35" - 48"	40" - 56"	50" - 65"	54" - 74"
Target User Age	1 - 2 yrs	2 - 6 yrs	6 - 9 yrs	9 - 14 yrs	13 yrs - Adult
Inside Seat Depth	7"	9.5"	11.5"	13.5"	16"
Inside Seat Width	8"	10"	12"	14"	16"
Inside Seat Height	21"	26"	29"	35"	39"
Floor To Seat Ht. Upright on Stationary Base	10"	10"	9.75"	11.5"	12.5"
Floor To Seat Ht. Tilted on Stationary Base	10.5"	10.5"	10.5"	13"	14.25"
Floor To Seat Ht. Upright on Mobile Base	12.5"	13.25"	13"	14.75"	16"
Floor To Seat Ht. Tilted on Mobile Base	12.75"	13.75"	13.75"	16"	17.75"
Floor To Seat Ht. Upright on Floor Base	5"	5.5"	5.5"	DNA	DNA
Floor To Seat Ht. Tilted on Floor Base	4.5"	4.5"	5"	DNA	DNA
Overall Seat Depth	12"	14"	16"	19"	21"
Overall Seat Depth	12"	14"	16"	18"	20"
Overall Seat Depth	25"	29"	34"	40"	44"
Weight of Sitter	6 lbs.	10 lbs.	13 lbs.	15 lbs.	19 lbs.