

# Out & About Seat (OAS)

Product Number: 67010100

MSRP: \$990.46

Provides comfortable support in a portable package

- Padded chest & pelvic harness hold your child securely in place allowing them to discover their world.
- Optional Special Tomato seat & back liner cushions enhance lateral support and comfort.
- Use in classrooms, restaurants, ball parks, movie theaters, relatives homes when visiting
- Assists your child for participation in classroom & social settings.
- Provides comfortable upper body support.
- Backpack straps keep your hands free making traveling a breeze.
- Even folds with Special tomato seat & back liners attached.
- Lightweight, portable and easy to clean

**Ki** Mobility  
Canada



*Standard OAS*





*Ready for travel as a backpack*



Toll-Free: (877) 252-6655 Fax: (403) 253-1744  
Email: [sales@kimobilitycanada.ca](mailto:sales@kimobilitycanada.ca)  
Order Forms Available Online at: [www.kimobilitycanada.ca](http://www.kimobilitycanada.ca)

## Out & About Seat Specifications


Special Tomato® OAS Size Chart	Special Tomato® OAS Dimensions
Seat Width	15 in
Seat Depth	13 in
Seat Back Height	27 in
Seat Back Height with Seat Liners	25 in
Minimum Chair Back Height	14 in
Seat base to top of head support Max	25 in
Shoulder Harness Height Levels	11.5 in, 14 in, 16.5 in, 19 in, 21.5 in
Overall Dimensions (width x length x height)	16 in x 49 in (flat)
Folded Dimensions (width x length x height)	16 in x 27 in x 6 in
Maximum Weight Capacity	90 lbs
Weight of OAS	7 lbs



*The Out & About Seat accepts seat & back liners*

Seat Cushions	Cherry	Lilac	Aqua	Dark Gray	Inner Dimensions	Outer Dimensions	MSRP
Size 2 (LS2)	72040201	72040202	72040204	72040209	10"W x 10"D x 4"H	12"W x 13"D	\$ 400.31
Size 3 (LS3)	72040301	72040302	72040304	72040309	12"W x 12"D x 4.5"H	14"W x 15"D	\$ 457.18
Size 4 (LS4)	72040401	72040402	72040404	72040409	14" W x 14" D X 4.5" H	17.5" W x 17.5" D	\$ 512.04

Back Cushions	Cherry	Lilac	Aqua	Dark Gray	Inner Dimensions	Outer Dimensions	MSRP
Size 1 (LB1) *	72050101	72050102	72050104	72050109	7.5"W x 10"H	10.5"W x 10"H	\$ 368.31
Size 2 (LB2) *	72050201	72050202	72050204	72050209	8.5"W x 13.5"H	12.25"W x 13.5"H	\$ 400.31
Size 3 (LB3) *	72050301	72050302	72050304	72050309	10"W x 14.5"H	13.25"W x 14.5"H	\$ 457.18
Size 4 (LB4) *	72050401	72050402	72050404	72050409	13" W x 15.5" H	14.5"W x 15.5" H	\$ 512.04



*\*Depth of Back Cushions are 1" thick at the narrowest point and 4" thick at the widest point*



Toll-Free: (877) 252-6655 Fax: (403) 253-1744  
 Email: sales@kimobilitycanada.ca  
 Order Forms Available Online at: www.kimobilitycanada.ca