

Please mark clearly with an X the items you wish to purchase.
 Email form to Ki Mobility Canada Ltd. at ca-sales@kimobility.com.
 Contact our customer service team at (877) 252-6655.

Account #: _____ Ship to: _____
 Date: _____ Name: _____
 PO #: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

Special Tomato Hi-Low™ Features:

- Affordable Customized Seating
- Tilt-in-Space with Child in Seat
- Floor to Table Adjustment
- Contoured Cushion Support
- Tool Free Adjustments
- Ergo Harness Design
- Soft-Touch® Material



Shown here with Standard Headrest (sold separately)

Adjustable 5-point harness

Shown here with Tray (sold separately)

Rear-locking casters

Adjustable Trunk cushions

Hip Flex™ Seat Cushion adjusts from 15° of hip flexion to 15° of hip extension



Select Your Cushions: Choose One of Each

MPS Seat Specifications	Small Seat Cushions				Large Seat Cushions			
	MS1	MS3	MS5	MS7	ML1	ML3	ML5	ML7
Target User Age	9 mos - 3 yrs	2 - 5 yrs	4 - 8 yrs	5 - 9 yrs	6 - 9 yrs	8 - 11 yrs	10 - 14 yrs	12 yrs - Adult
Target User Height	28" - 37"	32" - 43"	40" - 50"	43" - 52"	45" - 58"	48" - 60"	48" - 62"	58" - Adult
Maximum User Weight	80 lbs.				130 lbs.			
Seat Width	7"	9"	10"	10"	11"	12"	13"	13"
Seat Depth	6"	8"	10"	12"	11"	13"	15"	17"
Seat to Footrest	6" - 16"				8" - 18"			
Hip Flex Angle	-15° to +15°				-15° to +15°			
Seat Tilt Angle Forward	10°				10°			
Seat Tilt Angle Back	-30°				-30°			
Floor to Seat Height	14" - 28.5"				17" - 28.5"			

MPS Back Specifications	Small Back Cushions		Large Back Cushions	
	MS2	MS4	ML2	ML4
Target User Age	9 mos - 3 yrs	2 - 9 yrs	6 - 9 yrs	8 yrs - Adult
Target User Height	28" - 37"	32" - 52"	45" - 58"	48" - Adult
Maximum User Weight	80 lbs.		130 lbs.	
Back Width	8"	9"	10"	11"
Back Height	8.5"	11"	11"	13"
Top of Back Cushion	11" - 16"	15" - 20"	15" - 20"	19" - 22"
Top of Head Cushion	18" - 23"	21" - 26"	24" - 28"	26" - 30"

Tray Specifications	
Dimensions	24"W × 17"D × 2"H
Weight	8 lbs.

Toll-Free: (877) 252-6655 Email: ca-sales@kimobility.com
Order Forms Available Online at: www.kimobilitycanada.ca

Please mark clearly with an X the items you wish to purchase.
Email form to Ki Mobility Canada Ltd. at sales@kimobilitycanada.ca or fax form to 403-253-1744. Contact our customer service team at (877) 252-6655.

Customize Your Size



Small Shell (MS)

Shells
Choose a Shell size to determine your options



Large Shell (ML)

Options for (MS) Shell

Headrests

MHS (Standard)



7.5"W x 6.5"H


MHL (w/ Laterals)



7"W x 6"H

Back Cushions

MS2



8"W x 8.5"H

MS4



9"W x 11"H

ML2



10"W x 11"H


Seat Cushions

MS3



9"W x 8"D

MS5



10"W x 10"D

MS7



10"W x 12"D

Options for (ML) Shell

Headrests

MHS (Standard)



7.5"W x 6.5"H

MHL (w/ Laterals)



7"W x 6"H


Back Cushions

ML2



10"W x 11"H


ML4



11"W x 13"H


Seat Cushions

ML1




11"W x 11"D

ML3




12"W x 13"D

ML5



13"W x 15"D

ML7



13"W x 17"D

Please mark clearly with an X the items you wish to purchase.
 Email form to Ki Mobility Canada Ltd. at sales@kimobilitycanada.ca or fax form to 403-253-1744. Contact our customer service team at (877) 252-6655.



MPS Hi-Low™ Seating System

Includes Seat Shell with Footrest, 5-Point Harness, Hi-Low™ Base with Casters (Cushions Separate)

Size	Part #	MSRP
Small	79300109	\$ 5758.87
Large	79300309	\$ 7490.80



Optional MPS Tray

Sold Separately

Size	Part #	MSRP
Small	79060159	\$ 706.98
Large	79060359	\$ 797.29

Select MPS Cushions: One of Each



Headrest Cushions

Size	Cherry	Lilac	Aqua	Dark Gray	MSRP
MHS	79100051	79100052	79100054	79100059	\$ 259.80
MHL	79110051	79110052	79110054	79110059	\$ 323.72

Back Cushions

Size	Cherry	Lilac	Aqua	Dark Gray	MSRP
MS2	79120551	79120552	79120554	79120559	\$ 364.23
MS4	79120151	79120152	79120154	79120159	\$ 396.37
ML2	79120751	79120752	79120754	79120759	\$ 453.80
ML4	79120351	79120352	79120354	79120359	\$ 465.98

Seat Cushions

Size	Cherry	Lilac	Aqua	Dark Gray	MSRP
MS1	79050551	79050552	79050554	79050559	\$ 364.23
MS3	79050151	79050152	79050154	79050159	\$ 396.37
MS5	79050251	79050252	79050254	79050259	\$ 428.51
MS7	79050651	79050652	79050654	79050659	\$ 453.80
ML1	79050751	79050752	79050754	79050759	\$ 453.80
ML3	79050351	79050352	79050354	79050359	\$ 465.98
ML5	79050451	79050452	79050454	79050459	\$ 478.15
ML7	79050851	79050852	79050854	79050859	\$ 506.37

Does the intended user of this Hi-Low system weigh over 95 lbs. (43kg)? Yes No

\$

* Warranty:

5-Year Manufacturer's Warranty on cushions and seat shell
 2-Year Manufacturer's Warranty on Hi-Low™ base and strapping

Your MSRP Total

Toll-Free: (877) 252-6655 Email: ca-sales@kimobility.com
Order Forms Available Online at: www.kimobilitycanada.ca