

In-Situ Recline Deluxe

- Suitable for amputee scenarios
 - Based on the comfort of the Recline sling but with the added benefit of the Deluxe Leg sling design
 - Offers full body hammock style support
 - Includes additional leg paddles which provide femur support
 - Easy to apply from lying position with low friction fabric on leg sections to aid application
 - No aperture for optimum user comfort
- Soft binding as standard
 - Fixed colour coded Silva-Safety Slot® tapes
 - Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
 - Can be ordered with additional waist or chest support, head support and hip tapes if required
 - Wash temp 71°C / 159°F
 - Tumble dry low heat (max 50°C)

F15DLX-SFP



Product Codes & Specifications:



Silva-Superfine PLUS®		A	B	C	D
XS	F15DLX-SFP-XSM	76 - 86cm	81 - 92cm	72cm	97cm
S	F15DLX-SFP-SM	86 - 96cm	92 - 104cm	82cm	104cm
M	F15DLX-SFP-MED	96 - 108cm	104 - 116cm	92cm	110cm
L	F15DLX-SFP-LRG	108 - 130cm	116 - 128cm	102cm	117cm
XL	F15DLX-SFP-XLRG	130 - 144cm	128 - 140cm	112cm	127cm

