

In-Situ Recline HS

F16

- Suitable for amputee scenarios
 - Adjustable integral head support with removable stiffeners
 - Designed to provide closed leg full body support on transfer
 - Pommel strap for additional comfort and safety
 - To be applied in supine position
 - Loops supplied for provision for hip tapes
- Soft binding as standard
 - Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
 - Fixed colour-coded loop tapes as standard
 - Wash temp 71°C / 159°F
 - Tumble dry low heat (max 50°C)



NB: Silva-Safety Slot® tapes available upon request

Product Codes & Specifications:



Silva-Superfine PLUS®		A	B	C	D
XXS	F16SFPXXS	66 - 76cm	99 - 105cm	60cm	96cm
XS	F16SFPXS	76 - 86cm	105 - 112cm	70cm	112cm
S	F16SFPS	86 - 96cm	113 - 120cm	78cm	127cm
M	F16SFPM	96 - 108cm	120 - 128cm	88cm	138cm
L	F16SFPL	108 - 130cm	129 - 139cm	106cm	152cm
XL	F16SFPXL	130 - 144cm	138 - 150cm	117cm	162cm
XXL	F16SFPXXL	144 - 160cm	150 - 160cm	130cm	167cm

