





Silva Jolly Jumper

- Good supportive harness
- Only two hoisting straps attached to the top of the head support
- This product is supplied with adjustable sliders as standard
- Provides support to improve standing tolerance and independent transfers
- Used for rehabilitation purposes, therapy and inclusion in school and family activities
- Can be used on an overhead hoist or a mobile hoist
- Wash temp 85°C / 185°F
- Tumble dry low heat (max 50°C)

Product Codes & Specifications:



	100% Polyester	A	B	C	D
XXXS	SILJJXXXS	38-58cm	50 - 55cm	24cm	47cm
XXS	SILJJXXS	59-77cm	56 - 62cm	27cm	77cm
XS	SILJJXS	77-87cm	60 - 66cm	33cm	82cm
S	SILJJS	87-107cm	66 - 74cm	38cm	87cm
M	SILJJM	97-117cm	74 - 80cm	43cm	91cm
L	SILJJL	112-132cm	80 - 88cm	45cm	91m



SILJJ

