





Silva Balance Sling

- Fully supportive sling with full head support
- Manufactured with counterbalance straps to enable active patients more freedom of movement
- Fixed colour coded loop tapes as standard to the outer leg tapes
- Divided legs
- Waist support fitted as standard
- Manufactured in Silva-Superfine Plus® with soft binding as standard for greater user comfort
- Wash temp 71°C / 159°F
- Tumble dry low heat (max 50°C)

Product Codes & Specifications:



	Silva-Superfine PLUS®	A	B	C	D
XXXS	SILBF-XXXS	56 - 66cm	50 - 55cm	55cm	67cm
XXS	SILBFXXS	66 - 76cm	56 - 62cm	60cm	72cm
XS	SILBFXS	76 - 86cm	60 - 66cm	68cm	77cm
S	SILBFS	86 - 96cm	66 - 74cm	74cm	85cm
M	SILBFM	96 - 108cm	74 - 80cm	77cm	93cm
L	SILBF-L	108 - 130cm	80 - 88cm	80cm	97cm
XL	SILBF-XL	130 - 144cm	88 - 100cm	99cm	111cm



NB: Drawing shows Silva-Safety Slot® tapes, this product is now supplied with fixed colour coded loop tapes as standard.

SILBF



NB: Image shows Silva-Safety Slot® tapes, this product is now supplied with fixed colour coded loop tapes as standard.